

Denture Soft Liners

Your denture or partial denture has been relined with Lynal, a soft liner.

The liner will be quite soft for the first day or two and will then harden up with time. As it does, you may experience some rubbing of the gum tissue resulting in sore spots. If this occurs, please let us know so we may adjust your denture.

Cleaning:

Do not immerse the liner in bubbling cleaner solutions such as Efferdent or in bleach. To clean the liner, use a soft tooth brush with toothpaste or soap and water.

Replacement:

The lifespan of the liner seems to vary from patient to patient from one to 3 months. Keeping the liner clean will prolong its life before needing replacement.